




Inami & Co., Ltd.

OPERATION MICROSCOPES

GENERAL OPERATION MICROSCOPES





X-Y TRANSLATOR WITH RE-CENTERING

Horizontal fine adjustment of the microscope 45mm in all direction is possible with the motorized X-Y Translator with automatic re-centering.



ONE TOUCH BULB CHANGING SYSTEM

Our new one touch bulb replacement system permits the change of the bulbs without any inconvenient by sliding the fiber cable port to the other side easily.



AUTOMATIC ZOOM SYSTEM RATIO 1:6

According to the demands of the surgeon, we have finally developed an automatic zoom system. From 5X to 30X we supply a focused clear sharp image same as other operation microscopes.



FULL-FUNCTION FOOT SWITCH

Special designed Full-function footswitch allows the surgeon without hands over almost all functions, zoom magnifications, X-Y translator device, illumination switch etc.



**L-0980A
AUTOMATIC ZOOM OPERATION
MICROSCOPE WITH X-Y DEVICE**

2 STEPS VOLTAGE SELECTOR:

We have developed a new voltage selector system with 2 steps selector, so this will be available for unstable voltages.

- 100V - 120V AC
- 200V - 240V AC





SPECIAL COUPLING

BY using this special coupling, ENT, Plastic, Microendodontics or any other surgeries will be available.



L - 0980M MANUAL ZOOM OPERATION MICROSCOPE



UP/DOWN FOOT SWITCH



**L-0990
OPERATION MICROSCOPE**



**L-0970
MULTI-OPERATION MICROSCOPE**

SPECIAL COUPLING

BY using this special coupling, ENT, Plastic, Microendodontics or any other surgeries will be available.



**L-0940 SD
ZOOM PORTABLE MICROSCOPE**



General Purpose
By removing the 45° Inclined Adapter,
it can be used for ENT, plastic,
Microendodontics or any other surgeries.



**L-0548
45° INCLINED ADAPTER**

Convenient Portable Hard Case



K-1422 (OPTION)



L-0940SD + K-1422