



Inami & Co., Ltd.



Certified Management System
DIN EN ISO 9001
JIS Q 9001
EN ISO 13485
ISO 13485



INAMI BASIC &
REASONABLE INSTRUMENTS

IBRI

CATARACT SURGERY (ECCE & Phacoemulsification)

S
U
R
G
I
C
A
L
I
N
S
T
R
U
M
E
N
T
S



www.inami.co.jp/english

Still believe quality costs? Inami Basic & Reasonable Instruments

IBRI is our selection of world's standard ophthalmic surgical instruments series for international market. Our intention to launch this brand is to wipe off preconceived idea "quality costs."

We have specially developed a new product line only for this series with our strict guide line of quality control for MDs' total satisfaction.



Perfectly arrayed serration as a result of mastership



So precisely meshed that it makes no gap in between and holds a thread

FORCEPS



IS-310



IS-313



IS-310

Forceps, Tissue
Non-slip serrated handle for surgical gloves.

IS-313

Dressing 115mm



IS-311



IS-315



IS-311

Forceps, Tissue
Non-slip serrated handle for surgical gloves.

IS-315

Dressing 110mm



IS-312

Forceps, Tissue : slightly curved, 1x2 teeth
Non-slip serrated handle for surgical gloves.



S-331

Forceps, Corneal & splinter, THORPE



M-8RN



M-9RN



M-8RN

Forceps, Tying,
McPHERSON: straight, with extra
fine platform for Microsurgery

M-9RN

angled



M-177C

MICRO Colibri Forceps, with fine suturing
platform 0.12mm



M-11R
Forceps, IOL, SHEPARD



M-11RF
Forceps, IOL, SHEPARD, fine tip



IM-160R
IOL, Implant and Capsule Fragment Forceps GASKIN
Non-slip serrated handle for surgical gloves.



M-174
Forceps, for inserting folded IOL



M-1602
Forceps, IOL, folding, UNIVERSAL



M-1605
Forceps, for inserting folded IOL, MIYATA



M-167RA
Forceps, Circular Capsulorhexis, UTRATA



M-1610
Forceps, Circular Capsulorhexis, temporal approach



M-1611
Forceps, Circular Capsulorhexis , with round handle




M-1612
Inamura Small Incision Capsulorrhesis Forceps VI



MH-108



MH-109



MH-108
Haemostatic Forceps,
straight, serrated tip

MH-109
curved

SCISSORS



S-510SS



S-510SB



S-510CS



S-510CB



S-615S




S-615C



S-520S



S-520B



S-510SS
Scissors, Eye : straight, pointed

S-510SB
straight, blunt

S-510CS
curved, pointed

S-510CB
curved, blunt



S-615S
Scissors, Capsulotomy, VANNAS, straight


S-615C
curved




S-520S
Scissors, Tenotomy, WESTCOTT: slightly curved, pointed

S-520B
blunt

SPECULUM



M-2021N
Speculum: open ring model



M-135
Speculum, BARRAQUER, small size